

Personal Computer Control System

Higher performance type, ensuring easy operation and sophisticated data processing in Windows™ environments

- This system is an upgraded version of the stand-alone system.
It has the same test items but provides higher data processing capabilities of a wide variety.
- All the operations, including test parameter setting, test start and end, data acquisition, data processing, data saving, and reanalysis of stored data, are carried out by the personal computer.
- Functions of this system include processing various types of curves, overlaid presentation of more than one curve, and zooming.
- This system permits viscosity correction by the Bagley plot method and calculation of shear sensitivity (Dshear stress/Dshear rate), which provides information on the molecular weight distribution.



Shimadzu Doc

The personal computer and printer,
and the bench are not included in the standard setup.

Specifications

Testing Unit		Control Unit	
1 Test method: Constant-pressure extrusion system using weights.	7 Piston stroke: Effective stroke 15 mm Stroke detector Linear potentiometer Measuring precision $\pm 0.5\%$ (for 15 mm stroke) Resolution 0.01 mm	1 CPU: Single chip microcomputer, 16 bit	2 Memory size: 512 Kbyte ROM and 320 Kbyte RAM
2 Type of tests: Constant temperature test and constant heating-rate test	8 Maximum sample size: 1.5 cm ³	3 Display: 5.7" monochromatic LCD, touch panel type	4 Personal computer Interface: RS-232C
3 Extruding pressure range: CFT-5000 0.4903 MPa, 0.4903 MPa step selection CFT-1000 0.098 - 9.807 MPa 0.098 MPa step selection	9 Piston diameter: 11.282 mm (1 cm ² in cross sectional area)	5 Clock Incorporated	
Pressure precision Better than $\pm 1\%$ of the set value	10 Cylinder diameter: 11.329 mm		
4 Measuring range: CFT-5000 8×10^{-10} - 1×10^6 Pa · s CFT-1000 1.6×10^{-10} - 2×10^5 Pa · s *) for 0.5 ID x 1.0 mm length die **) for 1.0 ID x 1.0 mm length die	11 Standard dies (Interchangeable): 0.5 ID x 1.0 mm length 1.0 ID x 1.0 mm length 1.0 ID x 2.0 mm length 1.0 ID x 10.0 mm length		
5 Heating method: Electric, 500 W	12 Operational requirements: Power source 100 VAC, 50/60 Hz, 7 A max. Pneumatic source About 0.5 MPa, 15 Nl/min Temperature 10 - 35°C Humidity below 25°C in wet bulb reading 30 - 80% (No condensation)		
6 Test temperature: Range (Ambient + 20°C) - 400°C Heat sensor 100Ω platinum resistor Measuring precision error of heat sensor $\pm (0.3^\circ\text{C})$ Control precision $\pm 0.2^\circ\text{C}$ in stable state Heating rate 0.5 - 6.0°C/min (0.1°C step); (For constant heating-rate tests)	13 Approximate dimensions and weight: Rheometer main body 840 ^W x 650 ^D x 1,530 ^H mm, 150 kg Attachments 60 kg		

Stand-alone System Specifications		Personal Computer System Specifications
1 Parameters Selection of constant temperature method (hereafter called constant) or constant heating-rate method (hereafter called heating-rate), piston cross sectional area, piston drive condition, testing load, preheating period, sample name, sample number, test temperature (Constant), number of test repeats (Constant), density, calculation start/end positions (constant), test start temperature (Heating-rate), test final temperature (Heating-rate), heating rate (Heating-rate), interval of flow rate measurement (Heating-rate), die dimensions (orifice diameter and length), selection of flow rate calculation methods (limiting method or automatic method), full scale time, full scale stroke		
2 Dataprocessing items Shear stress τ (Pa), flow rate Q (cm ³ /s), shear rate $\dot{\gamma}$ (s ⁻¹), viscosity η (Pa · s), melt flow rate (MFR)(*), softening temperature (Constant), flow beginning temperature (Heating-rate), 1/2 method, offset temperature method (Heating-rate)		
3 Graph generation Stroke-time curve (Constant), stroke-temperature curve (Heating-rate)		Viscosity-shear rate curve (Constant) Shear stress-shear rate curve (Constant) Viscosity-temperature curve (Constant) Shear rate-temperature curve (Constant) Viscosity-time curve (Constant) Shear rate-time curve (Constant) Viscosity-temperature curve (Heating-rate) Shear rate-temperature curve (Heating-rate) Overlay of up to 10 curves Zooming display (Constant and Heating-rate) Viscosity correction by the Bagley plot method Shear sensitivity N

(*). The cylinder dimensions of CFT-500DV100D system is different from that of the melt indexer based on ASTM/JSI standards, therefore, this measurement results are not always equal to those of the melt indexer.

Personal Computer Hardware Requirements

SHIMADZU CFT-500D - Dies

<i>Die Diameter (mm)</i>	<i>Die Length (mm)</i>	<i>Comment</i>
0.5	1.0	used for non filled plastics
1.0	1.0	New condition
1.0	5.0	used for non filled plastics
1.0	10.0	used for non filled plastics
1.0	10.0	used for Glass filled plastics
2.095	8.0	ASTM standard for MFR

This document was created with Win2PDF available at <http://www.win2pdf.com>.
The unregistered version of Win2PDF is for evaluation or non-commercial use only.
This page will not be added after purchasing Win2PDF.