STERISCorporation STERIS Canada
STERIS Canada Corporation

Type de document/Document Type

INSTRUCTIONSDocument N^o.
920-514-036Modèle/Model
Reliance® 380PG/680PG
Pharmaceutical Grade Washer

Titre/Title

**Temperature Transmitter's Calibration with Allen
Bradley Control**

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Revu et approuvé par/Reviewed and Approved by:Nom/NameSignature

André Lemieux, Eng.

Special Applications Engineering

DISTRIBUTION: - Manufacturing and Qualification Documentation Binders 920-514-015PRÉPARÉ PAR/PREPARED BY:
Simon Villeneuve
Jean-Guy Roy
Software programmer

QAF-0504-02

Rev. 3

1998-12-10

NOTE: This instruction is to explain how to proceed to do this particular job. It is also considered as a record to prove that the job was indeed done.

SERIAL # :		ORDER # :	
CLIENT EQ. # :			

NOTE:

CHANGES AND/OR TECHNICAL MODIFICATIONS:

CONFIDENTIAL

1. MATERIAL NEEDED

(Matériel requis)

Loop Calibrator

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Last calibration date: _____

STERIS #: _____

Block Calibrator #1

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Last calibration date: _____

STERIS #: _____

Block Calibrator #2 (if used)

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Last calibration date: _____

STERIS #: _____

2. IDENTIFICATION

(Identification)

Remove TEMPERATURE TRANSMITTER from washer and put it in a convenient place to perform calibration.

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Sump Temperature Transmitter (Reliance® 380PG/680PG)

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Serial number: _____

Output range: _____

Final Rinse Tank Temperature Transmitter (Option Reliance® 380PG/680PG)

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Serial number: _____

Output range: _____

Exhaust Temperature Transmitter (Option Reliance® 380PG/680PG)

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3. CALIBRATION VERIFICATION

(Vérification de la calibration)

Perform a 2 points accuracy verification for each transmitter. If read values are within acceptable range, calibration needs not to be performed.

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3.1 Verification

3.1.1 **Set Main disconnect switch to "ON" position.**

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3.1.2 **Wait until Operator Interface shows "LOGIN"**

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3.2.1 **Set temperature of block calibrator to 100 °F (37.7 °C)**

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If temperature difference (Delta) between Operator Interface and Block calibrator is outside the acceptable error on Low or High point, perform a calibration on this temperature transmitter.

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ACCURACY VERIFICATION MADE BY: _____ DATE: _____

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4.1 Setup

(Préparation)

4.1.1 Set Main disconnect switch to "OFF" position.

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4.1.2 Ensure that loop calibrator is in READ Mode and in 4-20mA Mode

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4.1.3 Connect loop calibrator at the proper location (DIAGRAM #1).

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4.1.4 Set Main disconnect switch to "ON" position.

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4.1.5 Wait until Operator Interface shows "LOGIN"

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
4.1.11 Press button of transmitter to be calibrated

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Transmitter	Model (Option)	Button
Sump temperature transmitter	Reliance® 380PG/680PG	SUMP
Final rinse tank temperature transmitter	Reliance® 380PG/680PG (Option)	TANK
Exhaust temperature transmitter	Reliance® 380PG/680PG (Option)	DRYING

DIAGRAM #1
Loop Calibrator wiring

Analog Input Card

Slot #6, I In 0 +:	Sump Temperature Transmitter	TB3-101
Slot #6, I In 1 +:	Final Rinse Tank Temperature Transmitter	TB3-102
Slot #6, I In 2 +:	Exhaust Temperature Transmitter	TB3-103

5. CALIBRATION

Fill table #2 with all read and modified values during procedure)

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5.1.1 **Set temperature of block calibrator to low set point: 100 °F (37.7 °C)**

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5.1.2 **Put TEMPERATURE TRANSMITTER inside block calibrator then let temperature stabilize on Operator Interface.**

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5.1.3 **Note temperature on Operator Interface (Table #2, A)**

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5.1.4 **Note Current in mA of loop calibrator (Table #2, B)**

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5.1.5 **On TEMPERATURE TRANSMITTER, turn "0" potentiometer to have x.xx mA on loop calibrator:**

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Transmitter	Model (Option)	mA
Sump temperature transmitter	Reliance® 380PG/680GP	9.33
Final rinse tank temperature transmitter	Reliance® 380PG/680PG (Option)	9.33
Exhaust temperature transmitter	Reliance® 380PG/680PG (Option)	9.33

5.1.6 **Press NUMERIC ENTRY button to enter low value. Enter 100°F (37.7°C) press ENTER then press DONE button.**

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5.1.7 **Note Temperature on Operator Interface (Table #2, C)**

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5.2.1 **Set temperature of block calibrator to 185 °F (85 °C) or put TEMPERATURE TRANSMITTER inside the second block calibrator with temperature already set at this value.**

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5.2.2 **With Temperature transmitter on block calibrator, Wait until temperature stabilizes on Operator Interface.**

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5.2.3 Note temperature of TEMPERATURE TRANSMITTER on Operator Interface (Table #2, D)

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5.2.4 Note current in mA of loop calibrator (Table #2, E)

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5.2.5 On TEMPERATURE TRANSMITTER, turn "SPAN" potentiometer to have xx.xx mA on loop calibrator:

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Transmitter	Model (Option)	mA
Sump temperature transmitter	Reliance® 380PG/680PG	13.86
Final rinse tank temperature transmitter	Reliance® 380PG/680PG (Option)	13.86
Exhaust temperature transmitter	Reliance® 380PG/680PG (Option)	13.86

5.2.6 Press NUMERIC ENTRY button to enter high value. Enter 185 °F (85 °C) press ENTER then press DONE button.

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5.2.7 Note Temperature on Operator Interface (Table #2, F)

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5.3 Accuracy Verification

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5.3.1 On Block Calibrator, set temperature to 142.3 °F (61.3°C)

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5.3.2 With Temperature transmitter on block calibrator, wait until temperature stabilizes on Operator Interface.

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5.3.3 Note Temperature on Operator Interface (Table #2, G)

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5.3.4 Note Temperature on block calibrator (Table #2, H)

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5.3.5 If temperature difference between Operator Interface and block calibrator is lower than 1.0 °F (0.5 °C), press DONE (table #2, I).

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If temperature difference between Operator Interface and block calibrator is greater than 1.0 °F (0.5 °C) contact the Service Technician (Table #2, I).

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5.3.6 Remove prints form printer and fold to this document.

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6. FORMULA 0 - 300°F (-17.8 – 148.9°C) TRANSMITTERS

6.1 Definition

Span : Maximum operating temperature - Minimum operating temperature
300°F (148.9°C) - 0°F (-17.8°C) = 300°F (166.7°C)
Low end of range: Minimum operating temperature = 0°F (-17.8°C)

6.2 Formula :

$$\text{mA output} = 16 \frac{(\text{Know value} - \text{low end of range})}{\text{Transmitter span}} + 4$$

6.3 Low Setpoint

Known value : **100°F (37.8°C)**
Low end of range : 0°F (-17.8°C)
Span : 300°F (166.7°C)

Formula :
$$\text{mA output} = 16 \frac{(100^\circ\text{F} - 0^\circ\text{F})}{300^\circ\text{F}} + 4 \quad \text{or} \quad 16 \frac{(37.8^\circ\text{C} - -17.8^\circ\text{C})}{166.7^\circ\text{C}} + 4$$

Results mA output: **9.33 mA**

6.4 High Setpoint

Known value : **185 °F (85 °C)**
Low End of range : 0°F (-17.8°C)
Span : 300°F (166.7°C)

Formula :
$$\text{mA output} = 16 \frac{(185^\circ\text{F} - 0^\circ\text{F})}{300^\circ\text{F}} + 4 \quad \text{or} \quad 16 \frac{(85^\circ\text{C} - -17.8^\circ\text{C})}{166.7^\circ\text{C}} + 4$$

Results mA output: **13.86 mA**